



Have you made the switch yet?

POSAinc® 360 Knee Adductor Testimonials

“With their **infinite adjustability** and the availability of multiple pad sizes, the POSAinc 360 Knee Adductors are ideal in allowing me to accommodate all of my clients, ranging from **pediatrics to large adults**. Adding to that, these adductors can be mounted to any manual or power wheelchair configuration with ease, which makes the POSAinc 360 Knee adductor **my clear number one choice.**”

Rich, CO (ATP)



“This product **saves me time and money** because I have less breakage than with other manufacturers that offer a swing-away option. The customers love them!”

John, TX (RTS)



“The POSAinc 360s are simple to use and go on easy! The assembly, installation, and adjustment took less than 15 minutes. The patients love them because they are so **user-friendly.**”

Iggy, IN (RTS)



“The hardware is well built and **handles the roughest of users** and remains in tact.”

Jon, WI (ATP, CRTS)

“The knee adductors have a tremendous amount of **flexibility** and **versatility**. I can have a quickly removable knee adductor mounted anywhere.”

David, MI (ATP)



“The POSAlinc 360 system **makes transfers much easier**. Even clients with limited upper extremity function can flip the adductors out of the way, offering greater independence. Finally, **a true innovation** in lower extremity positioning.”

Phil, WI (ATP)

“The infinite adjustment allows **precise pad placement** without having to replace bracketry.”
Mark, CA (ATP, CRTS)

“The therapists love that they stay on the chair and are not removable, so **patients can't lose them.**”
Morgan, MD (ATP)



“The 360 Swing-Away Knee Adductors have **solved many problems** in the seating field.”
Roger, MI (ATP)

“Once again, you guys have figured out how to make useful, **durable** gear that makes guys like me look good. Thanks.”
Vincent, PA (Adapted Equipment Specialist)

“It doesn't matter what chair type that I am using because the 360 can **swing away in any direction** so it is a versatile adductor.”
Jason, FL (ATP)

“One of the positioning issues I struggled with are legs that abduct out. I had tried thigh supports of all kinds, but what really was effective was pushing on the knee. It was difficult to find a product that extended out to the knee, was easy to attach and release for transfers. But then I tried the POSAlinc 360 Swing-Away Knee Adductors and it was **like magic!** They were **adjustable to any position** I was looking for and kept the knees together. Great product.”

Karen, MN (ATP)

POSAlinc 360 is sold exclusively by AEL (US Patent #8,573,549)

