



Permobil® Track Insert Kit

AEL now offers a Permobil Track Insert Kit to mount any AEL POSAlinc Knee Adductor to the Permobil UniTrack

HOW TO ORDER

1. Order a Track-Mount Knee Adductor

Choose between the POSAlinc 360 or POSAlinc Original

2. Order a Permobil Track Insert Kit (Part No. 11143)

Order 1 Kit per Knee Adductor

Permobil® Track Insert Kit

Part No. 11143

POSAlinc® 360 Knee Adductor

Swings Away in Any Direction

A twist of the release knob unlocks the mechanism, enabling the pad/post assembly to swing – in any direction

The patented pin/channel design ensures that the knee adductor is perfectly repositioned when the assembly has been reset

Infinite Adjustment

To maximize personal comfort and support, a wide variety of adjustments are possible with the signature POSAlinc post and link combination

Coding Friendly

Itemized part numbers are available for ease of coding under E0956+E1028

Lifetime Warranty

AEL stands by its hardware with a lifetime warranty

Free Shipping

All AEL POSAlinc 360 Knee Adductors ship US ground at no charge



POSAlinc® Original Knee Adductor

Infinite Adjustment

To maximize personal comfort and support, a wide variety of adjustments are possible with the signature POSAlinc post and link combination

Removable Option

Order a POSAlinc Original Knob Kit (as pictured) to provide your clients a tool-less solution for adductor removal and adjustment

Lifetime Warranty

AEL stands by its hardware with a lifetime warranty

Free Shipping

All AEL orders over \$100 ship US ground at no charge

